

Vigour Back Support

Lightweight yet tough to meet the daily lifestyle demands of active wheelchair users. Back support angle and seat depth offers exceptional adjustability and can be facilitated with the user in the wheelchair. Ventilated spacer-fabric and shell system increases air-flow through the back offering daily comfort.

VIGOUR LO

New invertible shell design means the Vigour Lo can be used either way up on the wheelchair. Sacral and lumbar height support allowing unrestricted upper body range of movement.

Designed for the individual who requires:

- A solid back support with the ability to self-propel
- Minimal postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- For clients that have mild torso asymmetries



VIGOUR MID

Lightweight and designed for efficient propelling, with scapular subtractions for unimpeded arm function.

Designed for the individual who requires:

- A solid back support with the ability to self-propel
- Minimal to moderate postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- In combination with lateral supports offers exceptional trunk and thoracic stability
- For clients that have mild torso asymmetries



VIGOUR HI

Solid full height back support available in three heights for maximum trunk and shoulder support.

Designed for the individual who requires:

- A solid back support
- Minimal to moderate postural support
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- In combination with lateral supports offers exceptional trunk and thoracic stability



MICRO-CUBE FOAM UPGRADE

Air-flow breathability for maximum control of temperature and moisture. Micro-cube cells provide zoned weight distribution enhancing comfort.

1209-2351-300





Dual-layer cover reduces **SHEAR FORCES** with air-mesh side panels to counter moisture and temperature

Add laterals as required from our **LARGE MODULAR RANGE**

Can be upgraded with **DIFFERENT FOAM PADDING OPTIONS** for excellent results

SLOTS & HOLE POSITIONS on shell allow easy installation of head supports, lateral supports and harnesses

SIZE CHART



- Red squares are standard sizes.
- Non-standard or custom sizes are available on request.

PRODUCT DETAILS

Width	Height	Capacity	Product Weight
10"	2"	200lbs (90kg)	Lo Back (16x10) 0.7kg
11"-15"	2.5"		Mid Back (16x16) 1.0kg
16"-17"	2.75"	300lbs (136kg)	Hi Back (16x16) 1.3kg
18"-20"	3"		

- Note: This is the weight of the back support only and does not include the hardware. For hardware weight, refer to last page.

PRODUCT PERFORMANCE

Postural Control	Red	Red	Red	4
Lateral Stability	Red	Red	3	
Postural Asymmetries	Red	Red	3	
Changing Needs	Red	Red	3	
Pressure Care	Red	Red	3	


- Note: This is a guide only and ratings are dependent on correct use of product.

- Vigour Lo Back Support
- Vigour Mid Back Support
- Vigour Hi Back Support

Codes for Standard Sized Back Supports

Our back support codes are meaningful!

Here's how to generate a product code to place your order.



EXAMPLE

A Spex® T Shape Back Support of width 18" and height 20", would have the code 1255-1820-OTS.

Stand out from the crowd!

Upgrade your cushion with a splash of colour, for no extra cost!

We offer 5 colour options for our Spex® back support outer covers alternative to the standard black fabric. Simply add the upgrade codes to your back support order.



Calypso Pink
1209-2491-028



Chilli Red
1209-2491-022



Granite Grey
1209-2491-090



Mango Orange
1209-2491-030



Ocean Blue
1209-2491-060


SPARE OUTER COVERS for any Spex® back support are available for purchase. Refer to the Spex® ordering and scripting guide for details, or get in touch with your local Spex® distributor.

Spex® Back Support Hardware Options

QUAD MOUNT (DEFAULT)


Quad Mount hardware is the default hardware option for Spex® Back Supports and uses four independent points of attachment to the back support canes of the wheelchair for ultimate strength and lateral positioning. This hardware type offers 40mm of seat depth adjustment, and is available in three different variants:

QUAD MOUNT (DEFAULT)




- Has a unique locking system that allows the back support to be removed easily
- Provides 10° of back support angle
- Weighs 0.9kg

FIXED MOUNT



- Back support is non-removable without the use of tools but is lighter (and a more cost-effective option)
- Limits adjustability and back support angle
- Weighs 0.6kg

HEAVY-DUTY QUICK-RELEASE




- Supports heavier users with stronger L-brackets
- Shares the benefits of quick release hardware
- Weighs 1.1kg

DUAL MOUNT (OPTION)

Dual Mount hardware has two points of attachment to the back support canes of a wheelchair and offers simpler set up, increased adjustability, angle, and growth. It utilises the same quick-release hook and latch operation to secure the back support. Please note, this hardware option cannot be used on most wheelchairs that have a horizontal back cross brace.

Dual Mount hardware allows for:


- 40mm of seat depth adjustment
- 30° of back support angle
- Weighs 1.1kg



SINGLE MOUNT (OPTION)

Single Mount hardware is a lightweight option for lower height back supports and has two points of attachment to the back support canes of a wheelchair. Like Dual Mount, this hardware option also offers simpler set up, increased adjustability, angle, and growth, and allows for:

- 45mm of seat depth adjustment
- 30° of back support angle
- Weighs 0.9kg




BACK SUPPORT HARDWARE ACCESSORIES/ATTACHMENTS

QUAD/SINGLE MOUNT WIDTH EXTENSIONS

Allows mounting of back supports up to 3" wider than the chair width, depending on the size of the back support. Cannot be used on Quad Mount HD.

- Weighs 0.1kg (set of 4)
- 1259-2490-100



QUAD MOUNT DEPTH EXTENSIONS

Increase back angle by 10° (set of 2 for upper mounts)

- Weighs 0.1kg
- 1259-2450-100

Increase seat depth by 4" (set of 4 for all mounts)

- Weighs 0.2kg
- 1259-2452-100

